



MoticEasyScan Pro 6

The MoticEasyScan Pro 6 is the best desktop slide scanner on the market today. Compact, powerful, and deployable anywhere, this scanner is the perfect assistant for any digital slide application. Share your slides with a single click, or switch to live mode for direct specimen observation. This scanning system sets a new industry benchmark with its affordable price and high performance, seamlessly integrating into any workflow and bringing lab-quality scanning to any work environment.

Objectives	Plan APOCHROMAT 20X/0.75 S Apo Objective 10X/0.3
Optical System	CCIS®
Scanning Time (15x15mm - full tissue)	Standard mode: 160s (40X - with 20X objective) Standard mode: 60s (20X - with 10X objective)
Resolution	40X: 0.26µm/pixel 20X: 0.52µm/pixel
Focusing Technique	Real-time autofocus
Scanning Camera	5.0 MP (2/3" high speed Sensor)
Nosepiece	3 hole
Light Source	10W LED (Lifetime: 25,000 Hours)
Slide Capacity	6 Slides (Single Tray)
Slide Tray	6 Slide Capacity
Slide Dimensions	76 x 26mm
Slide Tolerances (mm)	Length: +0/-1, Width: +0/-1
Scanning Mode	Normal (Real-time autofocus) High precision (High precision autofocus) EDF (Extended depth of field) Z-Stack (Three Dimensional stacking)
Barcode Support	1D: Interleaved 2 of 5, Code 39, Code 128 2D: Data Matrix, QR Code
Computer	Included: Minimum Specifications: Intel Core i7-7700 16GB Memory 128GB SSD & 1TB SATA Disk Windows 10 Professional 64-bit
Monitor	Included: All-in-One 23.8" LED 4k resolution
Interface	USB 3.0
Included Software	DSAssistant EasyScanner software (for MoticEasyScan One)
Optional Software	DSAConference, DSServer, FS-Live Telepathology System, MoticFlow
Optional Objective	S Apo Objective 40X/0.75
Optional Modules	Large Format Scanning
Dimensions	33 x 33 x 44cm
Net Weight	16kg